

Wave 2018

- 1. Measure** your client to determine the size you need.
(See dimension chart for help.)
- 2. Select** the color, appropriate bath chair and chest strap.
- 3. Choose** the accessories you need.
- 4. Add** up the total cost.
- 5. Order** online in minutes at **Rifton.com** or see the back of this catalog for other ways to order.
- 6. Call** us any time at **800.571.8198** if you need assistance.
We're always here to help.

The tool-free **conversion kit** converts your small chair to medium as your child grows.

Fabric without calf rest

Blue	Z268	\$310
Pink	Z269	\$310

(Includes seat, back, cover and chest strap)

Fabric with calf rest


Blue	Z266	\$385
Pink	Z267	\$385

(Includes seat, back, cover and chest strap)



		Z210 small	Z220 medium	Z230 large	
Required components	Chair	Z211 \$260 <input type="checkbox"/>	Z221 \$270 <input type="checkbox"/>	Z231 \$330 <input type="checkbox"/>	
	Fabric cover	CHOOSE Blue	Z244 \$185 <input type="checkbox"/>	Z246 \$195 <input type="checkbox"/>	Z248 \$205 <input type="checkbox"/>
		CHOOSE Pink	Z245 \$185 <input type="checkbox"/>	Z247 \$195 <input type="checkbox"/>	Z249 \$205 <input type="checkbox"/>
	Chest strap	CHOOSE With lateral positioning	Z212 \$62 <input type="checkbox"/>	Z222 \$72 <input type="checkbox"/>	Z232 \$82 <input type="checkbox"/>
CHOOSE Without lateral positioning		Z213 \$41 <input type="checkbox"/>	Z223 \$52 <input type="checkbox"/>	Z233 \$62 <input type="checkbox"/>	
Optional accessories	Head blocks	Z214 \$88 <input type="checkbox"/>	Z214 \$88 <input type="checkbox"/>	Z234 \$88 <input type="checkbox"/>	
	Leg straps	Z215 \$41 <input type="checkbox"/>	Z225 \$52 <input type="checkbox"/>	Z235 \$62 <input type="checkbox"/>	
	Calf rest (calf rest must be chosen at the time of the original purchase. It cannot be added later.)	Z216 \$115 <input type="checkbox"/>	Z226 \$130 <input type="checkbox"/>	Z236 \$145 <input type="checkbox"/>	
	Tub stand	Z217 \$260 <input type="checkbox"/>	Z217 \$260 <input type="checkbox"/>	Z217 \$260 <input type="checkbox"/>	
	Shower stand	Z218 \$700 <input type="checkbox"/>	Z218 \$700 <input type="checkbox"/>	Z218 \$700 <input type="checkbox"/>	
	Tub transfer base (calf rest recommended)	CHOOSE Counterclockwise into tub	Z227 \$825 <input type="checkbox"/>	Z227 \$825 <input type="checkbox"/>	Z227 \$825 <input type="checkbox"/>
		CHOOSE Clockwise into tub	Z228 \$825 <input type="checkbox"/>	Z228 \$825 <input type="checkbox"/>	Z228 \$825 <input type="checkbox"/>
	Legs for tub transfer base	CHOOSE Standard	Z224 \$52 <input type="checkbox"/>	Z224 \$52 <input type="checkbox"/>	Z224 \$52 <input type="checkbox"/>
		CHOOSE Long (tub edge over 20" high)	Z238 \$52 <input type="checkbox"/>	Z238 \$52 <input type="checkbox"/>	Z238 \$52 <input type="checkbox"/>
	Additional head blocks	Z214 \$88 <input type="checkbox"/>	Z214 \$88 <input type="checkbox"/>	Z234 \$88 <input type="checkbox"/>	
	Additional fabric cover	Z219 \$185 <input type="checkbox"/>	Z229 \$195 <input type="checkbox"/>	Z239 \$205 <input type="checkbox"/>	
	Additional chest strap	CHOOSE With lateral positioning	Z212 \$62 <input type="checkbox"/>	Z222 \$72 <input type="checkbox"/>	Z232 \$82 <input type="checkbox"/>
CHOOSE Without lateral positioning		Z213 \$41 <input type="checkbox"/>	Z223 \$52 <input type="checkbox"/>	Z233 \$62 <input type="checkbox"/>	
Total	Z210 \$ _____	Z220 \$ _____	Z230 \$ _____		

Dimensions of the Wave

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 Key user dimension: height Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth. Important: Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.			
Item dimensions (inches)	Z210 small	Z220 medium	Z230 large
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½–18½	12½–18½	12½–18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Width of base	18	18	21
Length of base when flat	29½	29½	29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lbs)	75	150	250
Back and seat angle adjustments		0°, 15°, 30°, 45°, 60°, 75°	



Without calf rest
(Z269)

With calf rest
(Z267)

The conversion kit (two styles shown above) is the smart solution for growing children. Convert the small to a medium bath chair for a fraction of the cost of a new chair.

We'll send you a new larger fabric cover, back and seat frames and chest strap with instructions on how to remove the old and install the new parts on your existing base.



The Wave has a washable fabric cover that is easily removed for cleaning.



All three sizes are available in blue and pink.